

Modular Treatment Cabinet with Two Drawers and Two Doors

Model 6214

The UMF Medical 6214 Modular Treatment Cabinet has two drawers, two doors, a stainless steel push handle and four carpet casters. It comes equipped with a Kydex work surface and is the perfect complement to any exam room or treatment area.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › 2" (5.1cm) twin wheel carpet casters
- › Full extension roller bearing drawer glides
- › Stainless steel push handle
- › Heavy duty concealed door hinges
- › 3-year limited warranty
- › Kydex top with rail

Specifications

Width: 25" (63.5cm)
Height: 34.25" (87cm)
Depth: 18" (45.7cm)
Net Weight: 88lb (39.1kg)

Shipping Weight: 102lb (46.3kg)
Shipping Dimensions:
L30xW19xH39 (76.2cm x 48.3cm x 50.8cm)

Accessories

122 – Locking swivel casters 2" (set of 4)
502 – Double door lock (factory installed)
505 – Single drawer lock (factory installed)
507 – Additional keys (x2)