

■ SELECT COLLECTION

The Care Cliner

Available: Manual, Manual XL



Our medical seating portfolio has been constructed to optimize the clinical journey for patients and their caregivers.

The care continuum extends across many different environments, so Champion designs products to serve a variety of needs.

With an all-inclusive design, The Care Cliner is an ideal workhorse for treatment within extended care, long-term care, and clinics.



Key Benefits

- 3-position recliner with Trendelenburg (attained by attendant only). Back release handle for quick Trendelenburg positioning. Handle can be operated from either side of the chair.
- Elastic support system in seat and back frame for enhanced comfort.
- Can be used in dialysis, phlebotomy, oncology, and cardiac settings.
- For quality engineering and patient-centric innovation that prioritizes functionality, safety, and comfort, Champion products comply with ISO 13485 and are made in the USA.

Common Applications:

Extended Care, CHCs, Residential Facilities

Standard Features



Trendelenburg Positioning



Stationary Arms & Headrest Protector




Dual Fold-Down Side Tables




5" (12.7 cm) Casters

Space Requirements


Front View



Top View



Side View



Specifications

	Standard	XL
Model Number	6530	6540
Seat Height	21" / 53.34 cm	22.50" / 57.15 cm
Seat Width	21" / 53.34 cm	25" / 63.50 cm
Seat Depth	20" / 50.80 cm	20" / 50.80 cm
Overall Height	49" / 124.46 cm	49" / 124.46 cm
Overall Width (tables down ♦ tables up)	32.25" ♦ 50.75" 81.92 cm ♦ 128.91 cm	36.25" ♦ 54.75" 92.08 cm ♦ 139.06
Back Height Above Seat	31.50" / 80.01 cm	31.50" / 80.01 cm
Overall Depth (upright)	39" / 99.06 cm	44" / 111.76 cm
Overall Depth (reclined)	76.50" / 194.31 cm	76.50" / 194.31 cm
Armrest Height	27" / 68.60 cm	27" / 68.60 cm
Product Weight	94 lbs. / 42 kg	111 lbs. / 50 kg
Weight Capacity	275 lbs. / 124 kg	450 lbs. / 204 kg


Available Options



Heat + Massage Interior Arm Switch



Heat System Interior Arm Switch



Foley Bag Hook Accessory



IV Pole Options

(Additional available options are also offered. Talk to a Sales Rep for details.)