

Overbed Table, Double Top

HCPCS:



FIGURE A



- Features two table surfaces - -a larger one that tilts in either direction (Figure A) and a smaller one that stays securely in place
- Easily assembled
- Walnut, wood grain low-pressure laminate top
- Tops can be raised or lowered in infinite settings between 30" - 46"
- Swivel casters allow for easy maneuverability
- Chrome-plated steel "H" base provides security and stability
- Tabletops lock securely when height adjustment handle is released
- Tabletop level can be raised with the light upward pressure

PRODUCTS	Description
13068BV	OVER BED TABLE, DOUBLE TOPS, BV, 1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	55 lbs
Base Dimensions	25" (W) x 28" (D)
Carton Shipping Weight	24.5 lbs
Height Adjustment	30" - 46"
Weight	22 lbs