



Minor Laceration Tray with Instruments

ITEM 4523 20 trays / cs

CurPro Minor Laceration Trays are convenient, pre-assembled kits equipped with the essential components to care for minor lacerations. These compact, disposable sterile kits are designed to help protect against cross-contamination.

- Convenient kits contain necessary components for routine care of minor lacerations
- Includes: 1) Plastic Tray; (5) Gauze Sponges; (1) Waterproof Drape, blue; (1) Fenestrated Waterproof Drape; (2) Absorbent Towels; (1) Tissue Forceps; (1) Needle Holder; (1) Scissors
- Kits are sterile and disposable to help protect against cross-infection
- Forceps are metal and are non-serrated
- Plastic tray provides a stable, sterile surface
- Removable Tyvek lid adds safety and convenience
- Sterile, not made with natural rubber latex

