







Dynarex Nitrile Gloves in a Bag are high-quality nitrile gloves designed for enhanced protection and tactile sensitivity. These ambidextrous, powder-free Gloves feature textured fingertips for increased grip and beaded cuffs to allow for easy donning.

2618	LG	1 pr/bg	500 bgs / cs
2619	XL	2 prs / bg	500 bgs / cs
2624	XL	1 pr/bg	500 bgs / cs
2626UB-2	ΧL	2 prs / bx	100 bxs / cs



- High-quality nitrile gloves crafted for optimal protection and comfort
- Designed for tactile sensitivity and durable puncture resistance
- Textured fingertips provide improved grip for enhanced performance
- Minimum thickness of 2 mils at the fingertip and palm
- · Secure fit around wrist and handconforming flexibility help prevent fatigue
- Beaded cuffs allow for easy donning
- Ambidextrous and powder-free
- Not made with natural rubber latex











