



PEN NEEDLES

Dynarex Pen Needles are sterile, hollow needles specifically designed to easily deliver injectable medications. The lubricated head maximizes patient comfort, even with daily use.

- Hollow needles designed to deliver injectable medications effectively and efficiently
- Lubricated to help improve patient comfort
- Compatible with most insulin pens, available in a variety of gauges and lengths
- Sterile, non-toxic, and non-pyrogenic
- Individually packaged for optimal hygiene
- Ideal for use in clinical and home healthcare settings
- Not made with natural rubber latex

7024	30G x 5mm	100 / bx	12 bxs / cs
7025	31G x 8mm	100 / bx	12 bxs / cs
7026	30G x 8mm	100 / bx	12 bxs / cs
7027	31G x 5mm	100 / bx	12 bxs / cs
7028	32G x 4mm	100 / bx	12 bxs / cs

