



VAGINAL SPECULUM

Disposable

Dynarex Vaginal Speculums are premium-quality disposable speculums designed to help provide safe, effective visualization during pelvic exams. These smooth, rounded Speculums are designed to adjust to a variety of positions for optimal convenience.

- Designed for safe, effective visualization during pelvic exams
- High-strength plastic construction provides optimal durability
- Adjusts and holds in a variety of positions for convenient one-hand operation
- Smooth, rounded design allows for improved patient comfort
- Single-use and disposable for optimal hygiene
- 4914-4916: Light-source adaptable (corded)
- Not made with natural rubber latex or DEHP

4911	Small	10 / bg	10 bgs / cs
4912	Medium	10 / bg	10 bgs / cs
4913	Large	10 / bg	10 bgs / cs
4914	Small, with Light Option	25 / bg	4 bgs / cs
4915	Medium, with Light Option	25 / bg	4 bgs / cs
4916	Large, with Light Option	20 / bg	4 bgs / cs

