

HEAVY DUTY TRANSPORT STRETCHER

The Gendron Heavy Duty Transport Stretcher is designed for treatment and transport of larger patients. Proudly made in the USA, this Heavy Duty Stretcher offers a weight capacity of 1,000 lbs evenly distributed. Gendron stretchers are designed with features that provide effortless and reliable patient handling.

Features

- Hydraulic Fowler backrest adjusts from 0-80°
- Steel powder-coat painted construction
- Two-position vertical side rails
- Heavy duty 8" casters with steer and brake functions
- 4" high-density pad with nylon cover
- Four IV sockets — one at each corner of stretcher

Product Info



Specifications

	overall Width		overall Length		overall Height		Patient Surface				Side Rail				Clearance floor to base		Caster diameter		Fowler Angle	Shipping Weight ³		Weight Capacity ¹ EVENLY DISTRIBUTED	
	inch	cm†	inch	cm†	inch	cm†	Width	Length	Width	Length	Length	Height ²	Length	Height ²	inch	cm†	inch	cm†	degrees	lb	kg†	lb	kg†
G1190	40	101.60	80	203.20	29.50	74.93	39	99.06	79	200.66	57.50	146.05	10.50	26.67	10	25.40	8	20.32	0-80	350	158.76	1,000	453.59

¹ maximum weight capacity - evenly distributed - refer to user manual ² Height above patient surface ³ approximate measurement & options may alter weight specs [†] standard rounding method used



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

•1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

**graham
field**

GF2000234RevC25