

PATIENT TRANSPORT



GENERAL DUTY STRETCHER

The Gendron General Duty Transport Stretcher is designed for use in the hospital or clinic. Proudly made in the USA, the General Duty G1050 provides secure and comfortable transport throughout the facility. Gendron stretchers are designed with features that provide effortless and reliable patient handling.

Features

- One-function articulation
 - Raise and lower head section
- Carbon steel with carbon steel sheet metal construction
- Two-position fold-down-type side rails
- Heavy-duty 8" casters with side-mounted lever brakes equipped with Steer/Brake functions
- 2" high-density litter pad with low-shear fabric
- Four IV sockets — one at each corner of stretcher
- IV pole storage under stretcher top for optional telescoping IV pole item 10122

Product Info



Specifications

	overall Width		overall Length		overall Height		Patient Surface				Side Rail				Clearance floor to base		Caster diameter		Fowler Angle	Shipping Weight ⁵		Weight Capacity ¹ EVENLY DISTRIBUTED	
	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	degrees	lb	kg†	lb	kg†
G1050	28	71.12	77	195.58	29.50	74.93	23	58.42	75	190.50	57	144.78	11	27.94	10.50	26.67	8	20.32	0 - 75	210	95.25	350	158.76

¹ maximum weight capacity - evenly distributed - refer to user manual ² Height above patient surface ³ approximate measurement & options may alter weight specs ⁴ standard rounding method used



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

graham
field

GF2000023RevC25