

## OptimAL, Aluminum 294 Capacity Unit Dose Medication Cart, Standard Package #AL294BOX

OptimAL Line medication carts are lightweight due to the all-aluminum construction. These carts are available to accommodate all types of medication packaging and offer a complete line of accessories to meet your storage and security needs. Harloff carts are shipped fully assembled and ready to use.



### Construction:

- Painted aluminum cabinet and drawers

### Features:

- Dual column medication cart with 294 unit dose medication box capacity, single width narcotics drawer
- BEST® brand lock on cabinet with four keys
- Illinois lock on narcotics drawer with four keys
- Painted pull-out shelf
- Ball-bearing, full extension, self-closing drawer slides
- Plastic replaceable top with integrated push handles
- Protective corner bumpers
- 5" ball-bearing premium rubber casters, full swivel, 2 with toe locking brake and 1 directional
- Hard-baked powder coat finish
- Dividers are available to accommodate different medication packaging

### Left Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 6 ea - 3.25" drawers with internal dividers
- 1 ea – 9.75" multipurpose storage drawer with adjustable dividers

### Right Column Drawer Configuration:

- 1 ea – 4" multipurpose storage drawer
- 1 ea – 9.75" locking narcotics drawer
- 2 ea – 9.75" multipurpose storage drawer

### Assembled Cart Dimensions:

- 45.75" H x 43.75" W x 24" D  
(116.2 cm H x 111.13 cm W x 61 cm D)
- Writing surface: 44" H (111.8 cm)
- Cart weight: 162 pounds (73.6 kg)
- Approx. shipping weight: 205 pounds (93.2 kg)
- Approx. shipping dims: 57" H x 47" W x 24" D  
(144.8 cm H x 119.4 cm W x 68.6 cm D)

### Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

<https://www.harloff.com/finish-options/>

### Warranty:

- 12-Year Manufacturer's Warranty on shipments to US and Canada, 5-Year Manufacturer's Warranty on all other International shipments.
- 1-Year Limited Warranty on any bought-out parts