



# MPA1824E05+MD18-ANS

**ALUMINUM**

## A-SERIES MINI WIDTH CART w/Five Drawers and Electronic Pushbutton Lock

### Features:

- Extremely light weight cart with critical grade aluminum construction, ideal for carts that are moved consistently throughout the day
- Curved drawer pulls, durable and elegant
- Fully welded uni-body aluminum construction, highest quality in the industry
- Heat-sealed durable scratch resistant powder coat finish, inside and out
- Ordered accessories are easily attached to pre-threaded holes with a screwdriver
- Industry only fully curved protective bumper guards your cart and walls from damage
- Self-closing, steel ball bearing, full extension drawer slides for easy access and smooth feel
- Removable plastic top
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Fully customizable drawer configurations at no extra charge
- Unmatched selection of locking options. Ten options include multiple locks of different types on the same cart, alarm access locks and badge access lock systems. View online at <https://www.harloff.com/locking-options/>
- Multiple alternative storage options for items that won't fit inside a traditional drawer
- Cart is delivered assembled
- 24" of vertical drawer space

#### Drawer Configuration:

- Three 3" drawers (7.6 cm), one 6" drawer (15.2 cm) and one 9" drawer (22.9 cm)

**Warranty:** 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

All dimensions are available in the "MA-Dims-18wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.

**Engineered & Manufactured in the U.S.**



**Color shown with a Light Gray body and Teal drawers. More than 20 additional colors are available.**

### Optional MD18-ANS Accessory Package (shown) Includes:

- 3-Tier Back Rail System (MD18-RBR3)
- (2) Aluminum Side Rails (SIDERAIL12)
- 5 Compartment Tilt Bin with Rail Clips (TILTBIN5RC)

Enjoy top-quality innovative storage and organization with the sleek, modern appearance of A-Series carts. Lightweight aluminum with optimal curb appeal, the A-Series comes with curved drawer pulls, self-closing drawer slides and our expanded lock options for the perfect balance of form and function. A-Series carts offer you unlimited choices.



View similar products online at  
<https://www.harloff.com/product-category/anesthesia-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
[www.harloff.com](http://www.harloff.com) • Ph: (719) 637-0300**