



**MR6K**

## MR-CONDITIONAL ANESTHESIA CART with Six Drawers and a Key Lock

### Features:

- Tested and certified by independent testing authority to 3 Tesla level
- Meets American Society for Testing and Materials (ASTM) International, Designation: F2503-05 definition of "MR-Conditional"
- Designed and built for safe use in the strong magnetic fields present in Magnetic Resonance Imaging (MRI) suites
- Cabinet painted gloss white with MRI identification sticker applied
- Low-ferrous, ball-bearing, full extension drawer slides
- Protective corner bumpers
- Pull-out shelf
- Low-ferrous 5" sealed ball-bearing casters; two swivel – two with brake, one directional
- BEST® brand lock on cabinet, with 2 keys
- Removable plastic top with integrated push handles
- Cart is delivered fully assembled with pre-threaded accessory mounting holes
- MR-Conditional accessories available
- Latex free
- Durable, heat sealed powder coat finish



**Engineered & Manufactured in the U.S.**



**Color shown is White with Beige drawers. More than 20 additional colors are available.**

### Drawer Configuration:

- 23.9" W x 17" D (60.8 cm x 43.2 cm)
- One 4" drawer (10.2 cm), two 3.25" drawers (8.3 cm), two 6.5" drawers (16.5 cm) and one 9.75" drawer (24.8 cm)

### Carts Specially Made for Magnetic Resonance Imaging Suites

Harloff offers a line of carts that meet American Society for Testing and Materials (ASTM) International Designation of MR-Conditional for use in a 3 Tesla or less environment of the Magnetic Resonance Imaging Suite (testing date 1-13-2010). Whether you need it for Anesthesia or Emergency, Harloff has the right MR-Conditional cart for your needs.

#### Assembled Cart Dimensions:

- Writing surface: 43.5" (110.5cm)
- 46" H x 31.8" W x 23.75" D (116.8 cm H x 80.8 cm W x 60.3 cm D)
- Cart net weight: 108 lbs (49.5 kg)
- Approximate shipping weight: 146 lbs (66.4 kg)
- Approximate shipping dims: 55" H x 37" W x 26" D (139.7cm H x 94cm W x 66cm D)

**Warranty:** 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts



**View similar products online at**

<https://www.harloff.com/product-category/mri-safe-carts/>

**Toll Free: (800) 433-4064 • [medicalsales@harloff.com](mailto:medicalsales@harloff.com)  
www.harloff.com • Ph: (719) 637-0300**