

PEDIATRIC & YOUTH STRETCHER



THE 360 DEGREE PROTECTION SOLUTION

416 HORIZON® PEDIATRIC STRETCHERS

The Pediatric Stretcher is ideal for the smallest patients who need the most protection. For infants to six year olds, it provides the safety and security of 360° patient protection and access.

- 325 lb (147.42 kg) Weight Capacity – EVENLY DISTRIBUTED¹
- Adjustable backrest from 0° – 45° for patient comfort
- Dual Side Foot Controls for easy height adjustment and 10° Trendelenburg positioning
- Extra-high side rails provide three positions to aid in patient safety, access, and transfer
- Removable end rails allow easy patient access
- Brake and Steer Controls accessible at all four corners
- Built-in O2 Holder and Telescoping IV Pole Included

♦ See Chart for Specifications

 **416D0ST** - Pediatric Product Info



* shown with side and end rail pads accessory (H00011P00-xxx)²



462 HORIZON® SERIES YOUTH STRETCHERS

The Youth Stretcher is ideal for patients who are too large for a crib or pediatric stretcher but still require 360° protection. It is intended for easy transport, examination and care of smaller patients in healthcare facilities.

- 500 lb (226.80 kg) Weight Capacity – EVENLY DISTRIBUTED¹
- Dual Side Foot Controls for easy height adjustment and 12° Trendelenburg positioning
- Pneumatic-Assisted backrest from 0° – 90° to aid patient comfort and safer staff use
- 3-position vertical side rails to aid in patient safety and transfer
- Brake and Steer Controls accessible at all four corners
- Center Track Steering for superb mobility
- Built-in O2 Holder and Telescoping IV Pole Included
- Comfort Mattress Pad (optional upgrade)

♦ See Chart for Specifications

 **462EMPST** - Youth Product Info



Models: 416D0ST pediatric 462EMPST youth

	inch	cm [†]	inch	cm [†]
Overall Width:	29.75	75.57	30.00	76.20
Overall Length:	59.25	150.50	83.00	210.82
Overall Height: minimum	32.25	81.92	28.50	72.39
Overall Height: maximum	39.00	99.06	37.00	93.98
Patient Surface: width	25.00	63.50	25.00	63.50
Patient Surface: length	56.00	142.24	75.00	193.04
Side Rail Length:	56.00	142.24	73.00	185.42
Side Rail Height:	17.75	45.09	16.00	40.64
End Rail Length:	23.00	58.42	23.00	58.42
End Rail Height:	17.75	45.09	16.00	40.64
Spindle Gap: minimum	3.00	9.53	2.50	10.16
Spindle Gap: maximum	3.37	8.56	4.00	6.35
Mattress Pad: thickness	3.00	7.62	3.00	7.62
Caster Diameter:	8.00	20.32	8.00	20.32

	features	features
Caster Type:	Tente®	Tente®
Hydraulic Controls:	dual side foot pedals	dual side foot pedals
Backrest Control:	manual	pneumatic assist
Fowler Backrest: movement	degrees [°]	degrees [°]
Trendelenburg: standard	0 - 45	0 - 90
Trendelenburg: reverse	10	12
Weight Capacity: ¹ EVENLY DISTRIBUTED	lb kg [†]	lb kg [†]
Shipping Weight: ⁸	325 147.42	500 226.80
	332 150.59	383 173.73

[†] multiple fabrics options available - replace 'xxx' with color code [†] standard rounding method used
¹ maximum weight capacity - evenly distributed - refer to user manual ⁸ approximate measurement
 * options available at additional cost - accessories sold separately



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2014, GF Health Products, Inc. All Rights Reserved.
 Hausted, and Horizon are registered trademarks of GF Health Products, Inc.
 Tente is a registered trademark of TENTE International GmbH.

grahamfield

GF1400048 Rev F24