

PreserveTech Med-Aire Edge Alternating Pressure & Low Air Loss Mattress Replacement System



- Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and micro-climate control
- Designed to help prevent, treat and heal pressure injuries in the home or long-term care setting
- Control unit offers 10-minute cycles when alternating and can also be used in static mode (exclusive to #14360-B)
- Control dial to optimize patient comfort levels by weight (exclusive to #14360-B)
- Advanced digital control unit offers 4 adjustable APM cycle times (10/15/20/25), and adjustable comfort weight settings for optimal immersion & patient satisfaction (exclusive to #14360-P)
- 9 deluxe clinical and care giver features: auto-firm, seat inflate, micro processing feedback, lock-out and more (exclusive to #14360-P)
- Audio and Visual alarm alerts when air pressure is low
- Digital control unit features 3 diagnostic audio and visual alarms to alert for issues with air pressure, the system, or power (exclusive to #14360-P)
- System includes 20 individual air cells to maximize immersion and envelopment and redistribute pressure with pillow feature for head stability and comfort.
- 8" firm foam perimeter for safe ingress and egress (getting into and out of the bed), patient transfers and power outage protection.
- 3" foam base for power outage protection
- 5" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable and machine washable quilted 4-way stretch PreserveTech cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Cover is treated with an antimicrobial to guard against degradation from microorganisms on the cover. The durable cover also tolerates cleaning with disinfectants, refer to user manual for suggested disinfectants.
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design
- 45db decibel level

PRODUCTS	Description
14360-B	Med Aire Edge Mattress w/basic pump kit
14360-P	Med Aire Edge Mattress w/prem. pump kit
SPECIFICATIONS	Description
Weight Capacity	450 lbs
Height	8" (mattress only 5" air cells on top of 3" foam base)
Length	80"
Weight	35 lbs (mattress only)
ACCESSORIES	Description
14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA