

## Quick 'N Easy Comfort Mattress

HCPCS: E0272-1,E0272-1



FIGURE A




FIGURE B



FIGURE C

- Durable high quality latex free polyurethane foam mattress for utmost comfort and support
- Non-slip bottom and 4 corner straps keep the mattress securely in place
- Two side carry handles make lifting and placement easier
- Compression technology creates a 14" x 14" x 37.5" box footprint creating convenience for transport and delivery
- Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15076	Quick N Easy Comfort Mattress,1EA
15076C	Cover,15076,1ea
SPECIFICATIONS	Description
Warranty	1 Year Limited
Weight Capacity	275 lbs
FireRating	CFR 16 Part 1633 Federal Fire Code
Cover Material	Fluid-Resistant Nylon
Dimensions	36" (W) x 80" (L) x 6" (H)
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]
PARTS	Description
15076C	Cover,15076,1ea, EA