Ventyv®







Ventyv®

Nitrile PF Exam

PLUS 5.0 Blue



Ventyv° Nitrile Powder-Free Plus 5.0 Blue (Kangaroo) Exam Glove provides outstanding barrier protection with the fit and feel of a traditional nitrile glove.



Kangaroo

















Reorder #	Size	Gloves/ Box ¹	Boxes/ Case
10333103	XS	200	10
10334103	S	200	10
10335103	М	200	10
10336103	L	200	10
10337103	XL	180	10
10338103	XXL	180	10

PRODUCT ATTRIBUTES

- Blue
- Nitrile²
- Powder-free
- Tested for Use with Chemotherapy Drugs⁴
- Fingertip Textured

PHYSICAL ATTRIBUTES³

- Length: 220 mm (XS, S)

230 mm (M, L, XL)

- Width: 70 ± 10 mm (XS)

80 ± 10 mm (S)

95 ± 10 mm (M)

110 ± 10 mm (L)

120 ± 10 mm (XL)

THICKNESS (Single Wall)

- Finger: 0.14 mm / 5 mils

- Palm: 0.10 mm / 3 mils

PHYSICAL PROPERTIES

- Tensile strength before aging: 48 MPa
- Tensile strength after aging: 51 MPa
- Elongation before aging: 612%
- Elongation after aging: 541%







Outer cartons made with 90% recycled material



2. This product is not made with Natural Rubber Latex

3. Physical Attributes meet or exceed ASTM standards for exam gloves

4. CAUTION: Gloves used for protection against chemotherapy drug exposure must be selected specifically for the type of chemicals used. Review material safety data sheets for the drug being used to determine the required level of protection. Contact Ventyv for specific testing information.



