Ventyv®

Nitrile PF Exam

PLUS 3.5 Black







Ventyv®

Nitrile PF Exam

PLUS 3.5 Black



Ventyv® Nitrile Powder-Free Plus 3.5 Black (Stallion) Exam Glove has superior barrier protection, offering resistance to Fentanyl and a wide variety of chemotherapy drugs. It is ideal for use in a variety of first responder, medical, and emergency services settings.



Stallion

Reorder #	Size	Gloves/ Box ¹	Boxes/ Case
10318107	XS	100	10
10319107	S	100	10
10320107	М	100	10
10321107	L	100	10
10322107	XL	100	10

PRODUCT ATTRIBUTES

- Black
- Nitrile²
- Powder-free
- Tested for Use with Chemotherapy Drugs⁴
- Tested for Use with Fentanyl⁴
- Textured

PHYSICAL ATTRIBUTES³

- Length: 220 mm (XS, S)

230 mm (M, L, XL)

- Width: 70 ± 10 mm (XS)

80 ± 10 mm (S)

95 ± 10 mm (M)

110 ± 10 mm (L)

120 ± 10 mm (XL)

THICKNESS (Single Wall)

- Finger: 0.09 mm / 4 mils

- Palm: 0.07 mm / 3 mils

PHYSICAL PROPERTIES

- Tensile strength before aging: 46 MPa
- Tensile strength after aging: 50 MPa
- Elongation before aging: 630%
- Elongation after aging: 559%







- 1 By weight
- 2. This product is not made with Natural Rubber Latex
- 3. Physical Attributes meet or exceed ASTM standards for exam gloves
- 4. CAUTION: Gloves used for protection against fentanyl and chemotherapy drug expsoure must be selected specifically for the type of chemicals used. Review material safety data sheets for the drug being used to determine the required level of protection. Contact Ventyv for specific testing information



