

Full Sized-50 Bin Max (MTK1-4B771/TFTFT)

 [waterloohealthcare.com/product/50-bin-max-mtk1-4b771tfft-2](https://www.waterloohealthcare.com/product/50-bin-max-mtk1-4b771tfft-2)



Waterloo Healthcare's full sized Best key locking medication cart with drawers and up to 50 individual patient bins. Two keys included for Best central lock. Shown with optional color drawer trim.

Right Drawers (Half Sized) from top down:

- 4" With separate Best narcotics lock (includes two keys)
- 7"

- 7"
- 10"

Left Drawers (Full Sized) from top down:

- Cassette #1 (Shown with 3" bins)
- Cassette #2 (Shown with 5" bins)
- Cassette #3 (Shown with 3" bins)
- Cassette #4 (Shown with 5" bins)
- Cassette #5 (Shown with 3" bins)

SKU: MTK1-4B771/TFTFT Category: Medication Carts- Patient Bin Tag: Medication Cart

Height (in)	43	Drawers (Qty)	4
Width (in)	39	Lock Type	Key
Depth (in)	23	Material	Heavy-Duty Steel
Weight (lbs)	187	Pull-out Shelf	Yes
Warranty	5 Years	Built-In Dividers	No
		Plastic Top	Yes
		Anti-Tip Bumper	No
		Push Handles	Ergonomic
		Caster Config	2 Brake, 1 Tracking, 1 Swivel
		Ball-Bearing Slides	Yes