



HOSPITAL PRIVACY SCREENS

Dynarex Hospital Privacy Screens are convenient, rollable screens designed to provide enhanced patient privacy during medical procedures and exams. Featuring antimicrobial-treated panels and a lightweight design, these portable screens are ideal for a variety of healthcare environments.

- Versatile hospital screens designed to provide quick and effective patient privacy
- Three-panel or four-panel designs offer customizable configurations for a range of healthcare settings
- Lightweight aluminum frame with 3-inch hooded casters for effortless transport
- Folds easily for convenient, compact storage
- Antimicrobial-treated material helps prevent the spread of germs
- Durable nylon panels are flame-retardant for optimal safety
- Easy-to-clean surfaces for enhanced hygienic maintenance

10505	3 Panels	1 / cs
10506	4 Panels	1 / cs



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |